## DESCRIPTION

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Genzyme Biosurgery  
A division of Genzyme Corporation  
500 Kendall Street  
Cambridge, MA 02142, USA

## Références

Adams, M.E. (1993). An analysis of clinical studies of the use of crosslinked hyaluronan, hylan, in the treatment of osteoarthritis. *J. Rheumatol.* 20 (Suppl. 39), 16-18.

Adams, M.E., Atkinson, M.H., Lussier, A., Schulz, J.I., Siminovitch, K.A., Wade, J.P. and Zummer, M. (1995). The role of viscosupplementation with hylan G-F 20 (Synvisc®) in the treatment of osteoarthritis of the knee. *Osteoarthritis and Cartilage* 3, 213-226.

Balazs, E.A. (1968). Viscoelastic properties of hyaluronic acid and biological lubrication. *Univ. Michigan Med. Ctr. Journal*, 255-259.

Balazs, E.A. (1969). Some aspects of the aging and radiation sensitivity of the intercellular matrix with special regard to hyaluronic acid in synovial fluid and vitreous. In *Thule International Symposium: Aging of Connective and Skeletal Tissue*, October 1-3, 1968, (Eds. Engle, A. and Larsson, T.), Nordiska Bokhandeln's Forlag, Stockholm 107-122.

Balazs, E.A. and Denlinger, J.L. (1989). Clinical uses of hyaluronan. In *The Biology of Hyaluronan* (Ciba Foundation Symposium #143) (Eds. Evered, D. and Whelan, J.), John Wiley & Sons, Chichester and New York, 265-280.

Balazs, E.A. and Denlinger, J.L. (1993). Viscosupplementation: A new concept in the treatment of osteoarthritis. *J. Rheumatol.* 20 (Suppl. 39), 3-9.

Balazs, E.A. and Leshchiner, E.A. (1989). Hyaluronan, its crosslinked derivative—hylan—and their medical applications. In *Cellulosics Utilization: Research and Rewards in Cellulosics* (Eds. Inagaki, H. and Phillips, G.O.), Elsevier Applied Science, New York, 233-241.

Balazs, E.A. and Larsen, N.E. (1997). Safety data from seven clinical studies. Genzyme Biosurgery Report BXR 9009.

Band, P., Goldman, A., Barbone, A., Reiner, K. and Balazs, E.A. (1995). Intraarticular distribution and residence time of hylan polymers. Materials Research Society, Spring Meeting, April 17-21, 1995, San Francisco, CA, 433 (abstract).

Conrozier T., Bertin P., Mathieu P., Charlot J., Bailleul F., Treves R., Vignon E., Chevalier X. (2003) Intra-articular injections of hylan G-F 20 in patients with symptomatic hip osteoarthritis: An open-label, multicentre, pilot study. *Clinical Experimental Rheumatology*, 21; 5; 605-610.

Larsen, N.E. and Goldman, A.I. (1997). Residence time of Synvisc® (hylan G-F 20) in the rabbit knee. Genzyme Biosurgery Report BXR 9010.

Larsen, N.E., Lombard, K.M., Parent, E.G. and Balazs, E.A. (1992). Effect of hylan on cartilage and chondrocyte cultures. *J. Orthop. Res.* 10, 23-32.

Lussier, A., Cividino, A.A., McFarlane, C.A., Olszynski, W.P., Potashner, W.J. and de Médicis, R. (1996). Viscosupplementation with hylan for the treatment of osteoarthritis: findings from clinical practice in Canada. *J. Rheumatol.* 23, 1579-1585.

Peyron, J. G. (1993). A new approach to the treatment of osteoarthritis: viscosupplementation. *Osteoarthritis and Cartilage* 1, 85-87.

Pozo, M.A., Balazs, E.A. and Belmonte, C. (1997). Reduction of sensory responses to passive movements of inflamed knee joints by hylan, a hyaluronan derivative. *Exp. Brain Res.* 116, 3-9.

Raynauld, J.P., Torrance, G.W., Band, P.A., Goldsmith, CH., Tugwell, P., Walker V., (2002). A prospective randomized, pragmatic, health outcomes trial evaluating the incorporation of hylan G-F 20 into the treatment paradigm for patients with knee osteoarthritis (Part 1 of 2): clinical results. *Osteoarthritis and Cartilage.* 10:506-517.

Scale, D., Wobig, M. and Wolpert, W. (1994). Viscosupplementation of osteo-arthritic knees with hylan: a treatment schedule study. *Curr. Ther. Res.* 55, 220-232.

Wobig, M., Dickhut, A., Maier, R. and Vetter, G. (1998). Viscosupplementation with hylan G-F 20: A 26-week controlled trial of efficacy and safety in the osteoarthritic knee. *Clin. Ther.* 20 (3), 410-423.

Wobig, M., Bach, G., Beks, P., Dickhut, A., Runzheimer, J., Schwieger, G., Vetter, G. and Balazs, E. (1999). The role of elastoviscosity in the efficacy of viscosupplementation for osteoarthritis of the knee: A comparison of hylan G-F 20 and a lower-molecular-weight hyaluronan. *Clin. Ther.* 21 (9), 1549-1562.

Wobig, M., Beks, P., Dickhut, A., Maier, R. and Vetter, G. (1999). Open-label multicenter trial of the safety and efficacy of viscosupplementation with hylan G-F 20 (Synvisc®) in primary osteoarthritis of the knee. *J. Clin. Rheumatol.* 5 (6), S24-S31.

Genzyme Biosurgery, A division of Genzyme Corporation

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© 2015 Genzyme Corporation
SYNVISC et GENZYME sont des marques déposées de Genzyme Corporation.
sanofi-aventis Canada Inc.