



TOUJEO®: ADMINISTERED WITH A PREFILLED PEN

TOUJEO® [insulin glargine 300 Units/mL] is available in the TOUJEO® SoloSTAR® and DoubleSTAR® disposable prefilled pens (450 Units/1.5 mL and 900 Units/3 mL, respectively).

To avoid dosing errors and potential overdose, TOUJEO® must never be drawn from the cartridge of the prefilled pen into a syringe



Standard insulin syringes are graduated for 100 Units/mL of insulin, but the TOUJEO® prefilled pens contain 300 Units/mL¹

The dose pointers of the TOUJEO® SoloSTAR® and DoubleSTAR® disposable prefilled pens show the number of units of TOUJEO® to be injected.

The TOUJEO® SoloSTAR® and DoubleSTAR® pens have been specifically designed for TOUJEO®, so no dose re-calculations are required.

- Prior to initiation of TOUJEO®, train patients on proper use and injection technique
- Training reduces the risk of administration errors such as needle sticks and incomplete dosing
- Remind patients to follow the “Instructions for Use” to correctly use the pen device and administer TOUJEO®

TOUJEO® is indicated for once-daily subcutaneous administration in the treatment of adult and pediatric patients (6 years of age and older), with diabetes mellitus who require basal (long-acting) insulin for glycemic control.

Please consult the Product Monograph at <http://products.sanofi.ca/en/toujeo-solostar.pdf> for contraindications, warnings, precautions, adverse reactions, interactions, dosing, and conditions of clinical use. The Product Monograph may also be obtained by calling 1.888.852.6887.

Reference: TOUJEO® SoloSTAR® and TOUJEO® DoubleSTAR® Product Monograph, sanofi-aventis Canada Inc., May 1, 2020.